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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

DESCRIPTION

The acupressure pillow is a practical tool for relaxation, pain relief and support for body regeneration. Equipped with numerous acupressure points, it works on the principle of traditional Chinese medicine, stimulating the flow of energy and improving blood circulation. Thanks to its ergonomic shape, it is ideal for massaging the neck, back, nape, legs and other parts of the body. High-quality materials ensure durability and comfort of use.

APPLICATION

- Relaxation and stress reduction
- Pain relief
- Improved circulation
- Rehabilitation support
- For athletes and active people
- Energy increase

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- 73 massage rosettes
- Diameter of the rosette: 3cm
- Number of spikes on the rosette: 27
- Length: 38cm
- Width: 14cm
- Height: 10cm

SAFETY INSTRUCTIONS

- The acupressure pillow is designed to massage and stimulate acupressure points. Do not use it in a way that may cause discomfort or injury.
- Do not apply too much pressure to the pillow, especially on sensitive areas of the body. Use it gently, adjusting the intensity to your needs.
- Do not use on damaged skin.
- Avoid using the pillow on open wounds, scratches, skin inflammations or other injuries.
- People with health problems such as skin conditions, circulatory disorders, chronic pain or pregnancy should consult a doctor before using the pillow.
- If you experience severe pain when using the pillow, stop using it immediately and consult a specialist.
- Do not use the cushion for too long in one session. The recommended time is 10 to 30 minutes, depending on the sensitivity of the body.
- Avoid contact with sensitive body parts
- The pillow should not be used on the eyes, face or other sensitive areas of the body.
- The product is not a toy. Children should only use the cushion under adult supervision to avoid accidental injury.
- When using a pillow, place it on a soft surface such as a bed or mat to increase comfort and reduce the risk of injury.
- Make sure that all parts of the cushion, including the acupressure points, are in good condition. If there is damage, such as cracks or deformation, stop using the product.
- The pillow should not be used near heat sources or on hot surfaces to prevent damage to the material.
- If you notice skin irritation or an allergic reaction, discontinue use of the pillow and consult your doctor.
- Do not use in situations that threaten stability.
- Do not use the pillow while driving, operating machinery or in other situations requiring full concentration and stability.

STORAGE, CLEANING

- The pillow should be stored in a dry and well-ventilated room to prevent moisture accumulation and mold growth.
- Do not leave the pillow in the sun for long periods of time as UV rays may damage the material or cause the colours to fade.
- Do not store the pillow near radiators, stoves or other heat sources that may affect its durability.
- If the pillow is not used for a long time, store it in a cover or plastic bag to protect it from dust and dirt.
- Store the pillow in its natural shape to avoid deforming the acupressure spikes or damaging the filling.
- Wipe the surface of the cushion with a damp cloth and mild detergent. Make sure not to use aggressive cleaning agents that may damage the material.
- If the pillow has a removable cover, wash it by hand or in the washing machine on a gentle cycle at a maximum temperature of 30°C.
- Do not immerse the entire pillow in water
- If the pillow has a permanent construction, avoid completely submerging it in water to avoid damaging the filling.
- After cleaning, leave the pillow to dry in a ventilated place, avoiding direct contact with a heater or sun. Make sure the pillow is completely dry before using it again.
- Clean your pillow regularly, especially if you use it daily, to prevent the build-up of dust, bacteria, and sweat.
- You can use a fabric disinfectant spray every now and then to keep your pillow hygienic. Make sure the product is safe for the material.
- Do not wash the pillow or cover in hot water as this may damage the structure of the acupressure spikes and the fabric.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

The product meets European and national requirements on the safety of devices and products.

We reserve the right to make changes to text, design and product data without notice.

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